

Eustachian Tube Balloon Dilation System

Navilloon-e, Catalog ver. 2



The way back to your normal life

 **Navilloon-e**

**Eustachian tube balloon dilation system
with the quality you can trust**

Easy and simple procedure at your finger tip

- Grab and Push
- Bend and Stretch
- No cut, No hassle

Navilloon-e

- The innovative system to cure Eustachian Tube Dysfunction (ETD).
- The intro guide could be bent upon the user's convenience.
- A wider rim at the opening of the intro guide enables to make a full inflation of the balloon.



Intro guide

The intro guide is bendable upon the user's convenience.



Balloon

The balloon inflates the eustachian tube safely.

Rotation button

The rotation button enables to change the angle of intro guide.



Inflation part

The inflation part is the inlet where the pressure is injected by the balloon inflator.

Balloon controller

The balloon controller moves the balloon forward or backward.




Features

No.	Total Length	Weight	Balloon diameter & length	Balloon dilation
MG-BC-0601E	282.9mm	31.2g	6Ø (after dilation) x 20mm	12atm (Max)

※ For size, there may be an error range of ± 10%.

How to prepare



Prepare the Balloon Inflator, Navilloon-e catheter, and Intro guide unsealed.




Fill the balloon inflator with saline solution.



Keep the inflator vertically. Remove the air left in the Balloon Inflator.



Connect the Balloon Inflator tube to the inflation inlet of the catheter body. (If necessary, connect the extension tube to the inflation inlet, and the catheter body).




Adjust the direction of the intro guide.


How to use




Place the edge of the intro guide at the entrance of eustachian tube with an assistance of endoscopy visual system.



Insert the balloon into the eustachian tube by pushing down the Balloon controller. Keep inserting the balloon until you see the black mark.



Start inflating the balloon by injecting pressure with the balloon inflator until 12 atm is reached. Hold the position for 2 minutes to make sure the eustachian tube is dilated completely.



After 2 minutes, deflate the balloon with the balloon inflator. After deflation, withdraw the balloon probe by pulling back the balloon controller.



Seoul Office 5th Fl., CC Bldg. 439, Bongeunsa-ro, Gangnam-gu, Seoul, Korea

Head Office 110-1 Donghwagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, Korea

Call Center T.+82-70-5130-7800 F. +82-2-3661-0120

www.megamedical.co.kr [overseas@megamedical.co.kr\(overseas\)](mailto:overseas@megamedical.co.kr)